


# Vision Granulometry

Granulometry, shape, grain index...

Macroscopic or microscopic analysis by image processing



## Vision Granulometry

Turnkey systems dedicated to material analysis in your development and production laboratories. Vision Granulometry systems are deployed in chemistry, metal, pharmacy, automotive, cosmetics, food processing industries:



- ▶ Granulometry by opening
- ▶ Non-contact electronic sieving
- ▶ Size and shape distribution
- ▶ Grain index
- ▶ Particulate contamination analysis

- ▶ Modular and upgradeable turnkey systems - from millimeter down to submicron scale
- ▶ Compliant with international industry standards: ASTM, ISO, NF, VDA, USP, JIS, DIN...
- ▶ Automated, customizable reports, data export
- ▶ Secure settings by access levels
- ▶ User friendly interface and easy to use
- ▶ Accurate, reliable and reproducible results

# Vision Granulometry












Turnkey systems for all your granulometric analyses

<b>Granix</b> Granulometry by opening - NF X 11-696 Non-contact electronic sieving	<b>Ellix</b> Characterization of object size distribution depending on their shape
--	--



Acquisition	<ul style="list-style-type: none"> <li>▶ Calibration function, background correction for brightfield illumination</li> <li>▶ Integrated acquisition equipment steering and setting, image enhancement</li> <li>▶ Analysis from direct acquisition (scanner, zoom, microscope...), or based on archived image, video or mapping (image mosaics) issued from external acquisition systems (SEM, optical microscope, scanner...)</li> </ul>	
Features	<ul style="list-style-type: none"> <li>▶ Creation of study templates</li> <li>▶ Measurements beyond the camera field of view</li> </ul>	<ul style="list-style-type: none"> <li>▶ Secured operating modes - 3 levels</li> <li>▶ Automatic or step-by-step exploration</li> </ul>
Detection methods	<ul style="list-style-type: none"> <li>▶ User-defined or predefined (plugin) detection methods               <ul style="list-style-type: none"> <li>• Thresholding adapted to the contrast of the images, to the size of the objects, to the homogeneity of the lighting</li> <li>• Artifacts removal through standard or user-defined filtering</li> </ul> </li> <li>▶ Field analysis processed in 2 seconds</li> </ul>	<ul style="list-style-type: none"> <li>▶ User-defined or predefined (plugin) detection methods               <ul style="list-style-type: none"> <li>• Thresholding adapted to the contrast of the images, to the size and color of the objects, to the homogeneity of the lighting</li> <li>• Artifacts removal through standard or user-defined filtering</li> <li>• Agglomerates unbundling</li> <li>• Size and shape filtering applied to detected objects</li> </ul> </li> <li>▶ Definition of categories based on geometrical criterion (length, area, perimeter...)</li> <li>▶ Restriction according to minimum-size and maximum-size criteria</li> </ul>
Measurement	<ul style="list-style-type: none"> <li>▶ Mean diameter, mode diameter, standard deviation, min diameter, max diameter, undersize, separation limit, D10, D25, D50, D75, D90</li> <li>▶ Standardized or user-defined diameter classes (standard sieves): ISO 565, NF X11-501, ASTM E11, DIN 4188, BS410, GOST 3584, UNI 2331, JIS Z8801</li> </ul>	<ul style="list-style-type: none"> <li>▶ Objects position (X, Y) and angle, major/minor axis, thickness, Feret<sub>min</sub>, Feret<sub>max</sub>, area, perimeter, equivalent diameter, orientation</li> <li>▶ Length, width, elongation of the models</li> <li>▶ Aspect ratio</li> <li>▶ Models: elliptical, rectangular, rhomboidal, circular, square</li> </ul>
Representations	<ul style="list-style-type: none"> <li>▶ Histogram, cumulative undersize and oversize curves, normal distribution, statistical and digital values</li> </ul>	<ul style="list-style-type: none"> <li>▶ Counting tables, distribution histogram, scatter plot diagram and directional rose, statistical values</li> </ul>
Results and reports	<ul style="list-style-type: none"> <li>▶ Creation of automatic and customized reports</li> <li>▶ Data export to spreadsheet software or simply by Copy/Paste; reports recorded in MVR and pdf format</li> <li>▶ Keep the images from each field and include them in the report</li> <li>▶ Global mappings (image mosaics) and individual images (field by field) export</li> </ul>	

➔ **To go further, discover our specific systems**






**Filtrex**  
**Particle & fiber counting and granulometry**  
 ▶ ISO 16232, VDA 19, USP 788, ISO 4406, NF L41-101 compliance

**Tamlab**  
**Granulometry of sieve aperture**  
 ▶ NF ISO3310-1, ISO 3310-2, ISO 9044 compliance

**GrainSizer**  
**Automated grain size measurement**  
 ▶ NF A04-102, ASTM E112 compliance

Vision Granulometry **turnkey systems** are offered as complete systems (**optics**, video and **processing**) but also integrates easily into both your existing acquisition equipment and your computing system.

1 $\mu\text{m}$	10 $\mu\text{m}$	1 mm	FILE
<p><b>ULTRA FINE</b></p> <ul style="list-style-type: none"> <li>▶ Counting from 2.5 <math>\mu\text{m}</math></li> <li>▶ Motorized XYZ microscope and camera</li> <li>▶ Motorized or encoded stage - option</li> <li>▶ Brightfield or darkfield illumination / reflected or transmitted light</li> </ul>	<p><b>SUPER FINE</b></p> <ul style="list-style-type: none"> <li>▶ Counting from 10 <math>\mu\text{m}</math></li> <li>▶ Zoom microscope and camera</li> <li>▶ Motorized or encoded stage - option</li> <li>▶ Brightfield or darkfield illumination / reflected or transmitted light</li> </ul>	<p><b>FINE</b></p> <ul style="list-style-type: none"> <li>▶ Counting from 100 <math>\mu\text{m}</math></li> <li>▶ High resolution A3/A4 scanner</li> <li>▶ Transmitted light</li> </ul>	<p><b>FILE</b></p> <ul style="list-style-type: none"> <li>▶ All resolutions - from millimeter down to submicron</li> <li>▶ Images issued by your own acquisition system (optical or electronic microscope, photos...)</li> </ul>
 <p>Sample and virtual mapping review (overviews, picture mosaics) with a joystick</p>	 <p>Compatible with Zeiss, Nikon, Olympus, Motic, Optika devices...</p>	<p><b>Native 64</b></p> <p>Pictures acquired and saved in <b>full resolution</b></p> <p>Future analysis or comparison without loss of data</p>	 <p><b>OS:</b> compatible with Windows 7, 8, 8.1 and 10</p> <p><b>RAM:</b> 8 to 16 GB</p> <p><b>Processor:</b> iCore 7 minimum</p> <p><b>Peripherals and interfaces:</b> USB3 port and/or PCI express slot</p>

## By choosing Microvision, you can rely on:

- ▶ Our expertise in object characterization on microscopic and macroscopic scales - over 25 years serving the industry and life sciences
- ▶ User friendly systems developed in partnership with the key industrial players
- ▶ Strong network of resellers and international tech support

## Support and services

- ▶ Maintenance contract, Technical assistance
- ▶ Advice and expertise, Training

