


EVERBEING

Model: EB-005

**Sub-Micron/RF
Micropositioner**

RF

EB-005



 **EVERBEING**
EVERBEING INT'L CORP.

EVERBEING

Model:EB-005

Sub-Micron/RF Micropositioner

Features

- Sub-Micron Internal Circuit Probing
- Linear Motion
- Choice Of Triaxial/Coaxial Tip Holder
- Four-Sided RF Probe Positioner
East/West/South/North
- DC to 40GHz~120GHz
- Multi Contact Wedge DC to 40GHz or 67GHz
- Calibration Substrate With Software
- 45 Degree Cable Connection
No Need For Elbow Connector
- 30 Degree E-Z Probing Inclination
- Repairable Probe Head
- Tungsten Tip Available
- Positioner with Cable Holder
- Backlash-Free Movement

Specification

- | | |
|--------------|---------------------------------|
| • X-Y-Z | 12mmx12mmx12mm
Linear Motion |
| • Lead Screw | 200TPI Thread/Inch |
| • Resolution | 0.35Micron |

Dimensions

- 90mmWx130mmDx90mmH
- Weight 1000 Grams

Accessories

- Triaxial Tip Holder
- Coaxial Tip Holder
- Magnet, Vacuum & Locking Screw Based
- Soft Probing Tilting Mechanism
- RF Theta Mechanism
- Semi-Rigid Cables
- Active Probes
- Ultrasonic Cutter
- Low Current/Capacitance Probes
- High Voltage Probes



5F, No. 341, SEC. 1, KUANGFU RD., HSINCHU, TAIWAN
Tel: 886-3-666-2000 Fax: 886-3-666-2200
Mobile: 886-9111-3-9111
probestation@yahoo.com
Skype: probestation

